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AGRICULTURE AND AGRIFOOD DEPARTMENT

MINISTRY OF PRIMARY RESOURCES AND TOURISM | BRUNEI DARUSSALAM

GUIDELINES FOR SAMPLE COLLECTION AND SUBMISSION

ANIMAL DIAGNOSTIC (ADG)

Veterinary Laboratory Services Unit
Division of Livestock Industry
Department of Agriculture and Agrifood

First Edition 2013
Second Edition 2017
Third Edition 2019
Fourth Edition 2020

Compiled by
Veterinary Laboratory Services Unit
Division of Livestock Industry
Department of Agriculture and Agrifood
Ministry of Primary Resources and Tourism
Brunei Darussalam

Copyright 2013

All rights reserved. No part of this book may be reproduced in any form or by any means without prior permission from the Veterinary Laboratory Services Unit, Division of Livestock Industry, Department of Agriculture and Agrifood, Ministry of Primary Resources and Tourism, Brunei Darussalam

Perpustakaan Dewan Bahasa dan Pustaka Brunei Cataloguing-in-publication

BRUNEI Darussalam, Veterinary Laboratory Services Unit
Guidelines for Sample Collection and Submission – Animal Diagnostic / Brunei Darussalam,
Veterinary Laboratory Services Unit - Bandar Seri Begawan: Jabatan Pertanian dan
Agrimakanan, 2013.

p. cm.

ISBN 978-99917-38-62-8 (Kulit Keras)

1. Veterinary diagnostics I. Title.

636.089 BRU (DDC22)

TABLE OF CONTENTS

FOREWORD	1
LIST OF ABBREVIATIONS	2
GENERAL GUIDELINES FOR SAMPLE COLLECTION AND SUBMISSION	3
PROCEDURES FOR SAMPEL COLLECTION FROM LIVE ANIMALS OR AT POST MORTERM	6
SAMPLE FROM LIVE ANIMALS	7
1. Blood	7
2. Faeces	8
3. Skin Scrapping	9
4. Genital and Urogenital Tracts (urine discharges and semen)	10
5. Eye (scrappings and tears)	10
6. Nasal and Mouth (nasopharyngeal swab)	11
7. Milk	11
HISTOLOGY / CYTOLOGY	12
1. Tissue	12
2. Swabs and Smears	12
GROSS PATHOLOGICAL EXAMINATION (POST-MORTEM)	13
PARASITOLOGY EXAMINATION	14
1. Fresh Faeces	14
2. Whole Blood	15
3. Skin Scrapping/ Hair	15
4. Soil	15
SEROLOGICAL EXAMINATION (ANTIBODY / ANTIGEN DETECTION)	16

VETERINARY MICROBIOLOGY	17
MOLECULAR BIOLOGY	21
CLINICAL CHEMISTRY/ HAEMATOLOGY	
23	
TABLE 1: GUIDELINES FOR TAKING SAMPLE SIZES OF ANIMAL POPULATION	
24	
SAMPLE RETENTION TIME	
26	
CLIENT CHARTER	26
NEW CLIENTS	27

FOREWORD

Alhamdulillah, this Guideline for Sample Collection and Submission was developed to provide farmers, importers/ exporters of livestock and veterinary products, with information on sample collection and submission, type of samples and quantity of samples required for specific laboratory testing under Animal Diagnostic Sub-unit in Veterinary Laboratory Services Unit (VLS), Division of Livestock Industry, under the Department of Agriculture and Agrifood.

The detailed procedures and requirements in this booklet are based on the ISO/IEC 17025 and test method procedures for diagnostic testing requirement. This booklet also contains the client charters of each testing provided by VLS for client reference.

We hope that the Guideline for Sample Collection and Submission Animal Diagnostic (ADG) booklet can be used as a guide for farmers, importers/ exporters and the public who wish to send sample for diagnostic testing in VLS.

Finally, I would like to applaud and congratulate the VLS unit for their great effort in publishing this guideline. May Allah Subhanahu wa Ta'ala bless you all.

Wassalam

A handwritten signature in blue ink, appearing to be 'H.M.A. Yuliah Binti Abdullah', written in a cursive style.

[HAJAH MASLIANA YULIAH BINTI ABDULLAH]

Acting Director of Agriculture and Agrifood

LIST OF ABBREVIATIONS

Abbreviation	Meaning
ADG	Animal Diagnostic
BPW	Buffered Peptone Water
DNA	Deoxyribonucleic Acid
EDTA	Ethylene diamine tetra-acetic acid
ELISA	Enzyme-Linked Immunosorbent Assay
HE	Haematoxylin and Eosin
PBS	Phosphate buffered saline
RODAC	Replicate Organism Detection and Counting
VLS	Veterinary Laboratory Services

GENERAL GUIDELINES FOR SAMPLE COLLECTION AND SUBMISSION

All samples should be sent to the VLS as soon as possible after collection.

Samples Delivery Requirement:

- a) All samples shall be delivered to VLS within the working hours stated below:

Mondays---Thursdays and Saturdays

8.00 am – 11.30 am*

1.30 pm – 3.00 pm*

**Not applicable during the month of Ramadhan*

- b) Any sample sent outside these working hours will not be entertained.
- c) Whenever a sample must be delivered outside these working hours, the sender shall make pre-arrangement with an authorized person in charge (PIC), at least one day prior to submission.
- d) The above criteria (a), (b) and (c) do not apply to any *Crisis samples*

All samples collected would be categorised into three types:

- a) **Survey samples:** These types of samples are collected to detect the presence of certain infection in an appeared non-diseased animal population. A sample size for most laboratory examinations should then be considered sufficient when an expected prevalence rate of 5% is used (please refer to Table 1 for setting sample size).
- b) **Diagnostic samples:** These types of samples are taken from either any reported or known to be diseased animal populations, of which up to a 50% of the animals might presumably be affected with the infection (please refer to Table 1 for sample size).

c) **Highly infectious samples:** These samples are only withdrawn from an animal population known or suspected to be affected with a highly infectious disease and/or highly pathogenic to humans. One single sample taken would be considered sufficient, however this shall be classified a special case, and sampling shall be done by authorised officers only.

These following guidelines should mostly be applicable for samples under categories (a) and (b).

All samples sent to the VLS should be in good conditions, neither should be having undergone autolysis nor showing significant physical changes difference from that when samples just taken.

Any other abnormalities, which might totally be unrelated to that as being suspected as causes of death or diseases, such as due to cross-contamination and/or having undergone autolysis, should always be avoided.

All samples of which might be perishable or might contain materials sensitive to extreme temperature should always be kept in a suitable container for protection. Senders should only use suitable containers e.g. a cool box and required accessories for packaging and storing samples. This is part of safety measure to prevent any leakages and cross-contaminations to happen during transporting them to the laboratory.

All samples should be properly packed; clearly and precisely labeled using permanent or water-proof markers.

A cool box used for packaging and storing samples during submission to the laboratory should be filled with ice packs or ice to maintain temperature not exceeding 8°C.

Senders should consult laboratory personnel prior to sample collection depending on examinations required. Most samples might need special cares during submission, such as they might need a specific transport medium and/ or refrigeration (Procedure for Sample Collection from live animals and at Post-Mortem).

As soon as you reach the laboratory, please register the samples at the sample receiving counter. Fill in a sample request form, which is readily provided at the counter. Make sure the form is completed with all appropriate and necessary information required.

All required information to fill in the **Animal Diagnostic Test Request Form (ADTRF)** shall be completed in **BLUE** ink **ONLY**. Incomplete ADTRF is not acceptable.

PROCEDURES FOR SAMPLE COLLECTION FROM LIVE ANIMALS OR AT POST-MORTEM

Tissues and/or parts of which are deriving from both live animals and carcasses may be collected both at farm and during post-mortem respectively, depending on laboratory examinations required. Most certainly appropriate cares to all samples should be practiced in accordance with types of laboratory tests employed, such as microbiological culture, parasitology, clinical chemistry, histopathology and immuno-histochemistry, and detection of proteins or genome of nucleic acid.

All sampling shall be taken either directly by veterinarians or by the trained veterinary health personnel under the supervision of the veterinarians.

Each piece of tissue should be placed in a fully secured and clearly labeled separate plastic bag or sterile screw-capped jar / bottle.

Each piece of tissue or swab taken for microbiological examinations (bacterial or viral identifications) should be taken as aseptically as possible and may be sent in a transport medium depending on examination required.

SAMPLES FROM LIVE ANIMALS

1. Blood

- Samples may be taken for clinical chemistry, haematology, microbiological culture, direct visualization for bacteria, viruses or protozoa (*e.g. Trypanosoma spp*), in which case it is usual to use anticoagulants, such as K_2 EDTA or heparin.
- Samples should be taken, as clearly as possible, by venipuncture from the jugular vein, brachial (wing) vein, mammary veins, auricularis or retro-orbitalis veins (*i.e.* small animals) or cardiac puncture.
- Sample may be taken by using syringe with needle or by needle with vacuum tube on a designated site where skin should be shaved (plucked) and rubbed with cotton swab soaked in 70% ethanol.
- When blood smear samples are required, both thick and thin blood smears should be prepared from fresh samples and either they are then directly air-dried or methanol-fixed.
- For any required tests conducted in serology and DNA technology laboratories, a few drops of blood might alternatively be collected using filter papers. These samples should however be kept in low temperatures not exceeding 8°C during submission. Consult with Veterinary Officer before proceeding.
- Any sample collected with anticoagulants, thorough mixing by only gentle agitation is necessary immediately after taking.

2. Faeces

- An amount of at least 10g freshly voided faeces should be collected. Samples taken for parasitology should be kept in a clean and/or sterile plastic bag and should be sent refrigerated (chilled) using frozen ice packs to prevent hatching of parasite eggs (ova) when the temperature reaching up to 23° C.
- For bacterial isolation particularly *Salmonella*, sample collection should be taken from either fresh dropping (caecal) or from swabbing cloaca. Cares to prevent any cross contamination to the samples, an aseptic method of collection should be practiced. Such a sample should be then placed in a sterile polycarbonate jar/bag and kept refrigerated (chilled) using frozen ice packs.
- When *Giardia spp.* is suspected, samples can be placed in polyvinyl alcohol fixative at a ratio of one part faeces to two parts of polyvinyl alcohol.

3. Skin Scraping

- For viral isolation, approximately 2 grams of affected epithelial tissue (vascular lesions if any) should be taken aseptically and placed in 5ml phosphate-buffered glycerine (pH 7.6).
- Plucked hair or wool samples may be collected as samples, of which are useful material for identification of surface-feeding mites, lice and fungal infection.
- Skin (deep skin) scrapings of mammals should be collected using the edge of a scalpel blade for identification of burrowing mites; whereas feather tips of birds (avian) might be taken for detection of viral antigen, particularly when Marek's disease is suspected.
- All samples should be placed in bottle or sealed plastic bags without any preservatives added.

4. Genital and Urogenital tracts (urines, discharges and semen)

- Samples deriving from vagina and prepuce should be collected by either washing/flushing or suitable swabbing.
- Samples from the cervix or urethra should only be taken by swabbing.
- Urine samples should be freshly voided or quickly collected when target animals naturally pass the urine. If sterile sample is required, a cystocentesis method may be undertaken by a qualified veterinarian.
- Semen samples should be obtained either by using an artificial vagina or an electro-ejaculator. The samples should be taken care aseptically, placed in sterile tubes and sent in suitable transport medium.

5. Eye (scrapings and tears)

- Conjunctival swabs should be taken by holding the palpebra apart and gently rubbing the surface with dacron, cotton or gauze swabs. The samples are then put into a suitable transport medium.
- A mucosa scraping may also be taken by using a metal-handled sterile swab and then smeared the scraping onto a microscope slide for microscopic examination.
- 'Mucopurulent lacrimal' discharges may be taken as samples but rarely useful.

6. Nasal and mouth (nasopharyngeal swab)

- Nasopharyngeal swab should be taken for use in tests when viral infection is suspected, particularly in large ruminants.
- The sample may be taken using dacron, cotton or gauze swabs with metal wire handles. Swabs are suggested to be firstly moistened with transport medium, and then be allowed to remain in contact with the secretions for at least one minute. The swabs should be sent immediately and be refrigerated (chilled) while transporting.

7. Milk

- Samples should only be collected after cleaning and drying tips of the tits. The use of antiseptic for cleaning should however be avoided.
- When sampling, the initial stream of milk should be discarded, as only the following streams should be collected for samples.
- Samples for use in serological tests should not be frozen, heated or subjected to violent shaking. Adding permissible preservative to the samples for such a purpose might however be allowed, especially when there is going to be a delay in submitting them to the laboratory.

HISTOLOGY / CYTOLOGY

1. Tissue

- Sample collected for further use in microbiological examinations, particularly for bacterial and/or viral isolation, should be taken as aseptic as possible, and then placed in a sterile container with or without transport medium. The sample should be kept refrigerated (chilled) and sent immediately.
- Sample for further use in histopathological examination, particularly when only employing routine staining method e.g. HE, should be taken in blocks of tissue not more than 0.5 cm thick and 1.0 to 2.0 cm long, fixed in 4-10% formal-saline and placed in suitable jars. When a cryostat technique is required in order to directly visualize the presence of pathogens or inclusion bodies left by viruses, tissue should be taken as aseptically as possible and should immediately be frozen at -20°C using dry ice.

2. Swabs and smears

- Samples for use in microbiological examinations may only be required by swabbing main affected areas of carcass, especially when bleeding from ears, nostrils and mouth are indicated in large ruminants, and anthrax is then suspected, post-mortem inspection shall not be performed until a rapid test on swab results (Gram Staining) is negative for *Bacillus anthracis* presence.
- Impression smears may also be taken for direct visualization of pathogens.
- Samples shall immediately be sent to the laboratory, and always be regarded as being highly pathogenic ones.

GROSS PATHOLOGICAL EXAMINATION (POST MORTEM)

The whole carcass, with size not larger than an adult male deer, shall be sent to the laboratory immediately after death, not exceeding six (6) hours. In the events of postponement or a delivery cannot be done immediately after death, any carcass should be kept chilled or frozen prior to submission in the next following working hours.

To identify certain disease in a population of animals raised intensively such as flocks of poultry, and when an expected proportion of positive is 50%, a sample size sent to the laboratory should be at least four (4) to five (5) animals, either sick or recently dead ones taken from the respective population at risk.

Each individual of samples as the same size as chickens should be submitted in separated sturdy boxes and labeled appropriately.

PARASITOLOGY EXAMINATION

Current parasitological examinations carried out are only for the diagnosis of internal parasitism by employing methods for the detection and enumeration of worm eggs present in fecal samples of goats, sheep, buffaloes and cattle.

Faecal samples must be directly taken from rectum of target animals either by means of manual rectal exploration or using a rectal-probe. As much of 10 gram faecal sample is required, individually put in clean plastic bag, and properly labeled. The number of samples should be sufficient as representatives for respective animal populations (*please refer to Table 1*).

All samples should be stored in a cool box with ice packs to control temperature not exceeding 8°C during transporting to the laboratory and reach the laboratory not later than 6 hours after the sample collection.

1. Fresh faeces

- The sample should be taken and sent to the laboratory as immediately as possible (never held any sample more than 48 hours).
- Each sample should be placed in individual sealed containers, labeled properly and identified with animal number/ name and date of collection.
- Each sample should be kept chilled (not frozen) as soon as after collection.
- Each sample should be collected and sent at least 5 grams in quantity.
- For larvae examination (*e.g.* Baerman method) samples must be less than 2 hours old and unrefrigerated.

2. Whole Blood

- Samples with anticoagulant should be collected and sent at least 1 mL blood for example in K₂ EDTA tube or others depending on tests requested.
- Serum obtained from whole blood collection taken in a red cap tube. At least 0.5 mL should be sent to the laboratory in refrigeration (chilled with frozen ice blocks)

3. Skin Scrapping/ Hair

- Sample should be kept in a glass tube immediately after collection (a red-top blood collecting tube can be used).
- In case of small samples size, samples can be submitted on a glass slide if the cover slip is securely sealed with “Nail Polish”. Please avoid using a scotch tape as a cover slip because it will be very difficult or even impossible to observe detail through the tape.

4. Soil

- Soil samples should be collected at least 30 g each.
- The samples should be collected from two to three inches below the surface from at least three sites in the suspected area.
- Please send each sample in each sealed plastic zip lock bag.

SEROLOGICAL EXAMINATION (ANTIBODY / ANTIGEN DETECTION)

For ELISA tests, whole blood samples and / or serum (deriving from the target blood animals) should be immediately transported to the lab after they are collected. The samples can either be put in vacutainer blood tubes (provided by Serology Laboratory upon request) or in syringes. The volume of each sample will depend on the size of the target animals; ranges from 3.0 mL to 10.0 mL. Meanwhile, the number of samples taken shall be sufficient to represent the whole animal population (or a cohort) in one area, as required by the standard procedure used in animal disease survey (refer to Table 1). This procedure also applies to when collecting cloacal swabs, body fluids or hatching eggs.

For rapid test such as for Avian Influenza screening, both disposable swabs and assay diluents will be provided by the laboratory. The diluent of swab samples should be mixed well. After collection, the diluents shall be kept in a cool box filled with ice packs before transporting them to the laboratory.

VETERINARY MICROBIOLOGY

Test to be performed	Appropriate samples	Collection of samples
Isolation and Identification of <i>Salmonella</i> species from the farm environment	<ul style="list-style-type: none"> • Sampling with sterile swab tips (Surface sampling including ventilation ducts, drains or feed containers) 	<ul style="list-style-type: none"> ➤ Remove a swab from the sterile wrapping and moisten the tip by immersing it in a tube containing BPW. Press the tip of the swab against the wall of the tube to remove excess fluid. ➤ When sampling surfaces, a sufficiently large area should be swabbed to ensure that all surfaces of the swab are liberally coated with material. Ideally several different locations should be swabbed or multiple swabs taken to maximize recovery of the target organism. After sampling, break or cut off the stick aseptically. Place in transport medium (bottle containing BPW).
	<ul style="list-style-type: none"> • Sampling with boot swabs (Sampling floors) 	<ul style="list-style-type: none"> ➤ Slip a disposable Plastic Boot Cover over each shoe or boot after enter the house. ➤ Remove the pre-moistened boot swab from the Whirl-Pak[®] bag and place it securely over the plastic boot cover. ➤ Reseal the bag immediately after removing the boot swab to prevent contamination from airborne debris or insects. ➤ Step at least 100 steps in per house, including the full length and width of the house and any littered or slatted areas and ensuring that areas of faecal accumulation or wet litter under drinkers are included. ➤ Immediately after collection, carefully remove the boot swabs and return it to its original bag. Use appropriate measures to ensure there is a good seal on the Whirl-Pak[®] bag and wipe down the outside to remove any residual debris. ➤ Submit samples to the laboratory promptly for analysis.

Test to be performed	Appropriate samples	Collection of samples
	<ul style="list-style-type: none"> Litter samples (Sampling floors) 	<ul style="list-style-type: none"> The best way to collect litter samples is to move over the whole house taking pinches of litter from at least 60 separate areas in the house to make up a final amount of about 2 kg. Placed samples in a sterile bag. This sample can be sent to the laboratory or can be thoroughly mixed and a test portion of at least 25 g taken and sent to the laboratory. Litter samples should be taken before any replenishment of bedding.
	<ul style="list-style-type: none"> Water samples 	<ul style="list-style-type: none"> Collect at least 100 ml of water in sterile polypropylene sample containers / bottles with leak-proof lids (kept chilled during transport and storage). Use insulated containers to assure proper maintenance of storage temperature. Ensure that sample bottles are not totally immersed in water from melted ice during transit or storage.
	<ul style="list-style-type: none"> Animal feed 	<ul style="list-style-type: none"> Submit at least 500 g of sample. Chilling during transport and storage is not required.
Isolation and Identification of <i>Salmonella</i> species from animals	<ul style="list-style-type: none"> Cloacal swabs 	<ul style="list-style-type: none"> Collect using a sterile swab moistened with BPW. Ensure faecal material is present on the swab after collection. Break or cut off the stick aseptically. Placed in in a bottle containing BPW. The bottle with the samples is then well secured with the cap, labeled, surface decontaminated and placed in a cool box.
	<ul style="list-style-type: none"> Faeces 	<ul style="list-style-type: none"> Collect 10 to 20 g of fresh faeces in sterile bag and chilled. Fecal samples must not be held longer than a few hours. If they are held, they must be refrigerated. Do not submit feces in a glove.

Test to be performed	Appropriate samples	Collection of samples
Isolation and Identification of <i>Salmonella</i> species from hatchery	<ul style="list-style-type: none"> Broken eggshells 	<ul style="list-style-type: none"> ➤ Broken shells from hatched eggs can be collected from systems where basket liners are not in used. Egg fragments from as many baskets as possible per hatcher should be collected in a large bag.
	<ul style="list-style-type: none"> Dead-in-shell embryos (infertile eggs) 	<ul style="list-style-type: none"> • A minimum of 60 eggs is taken and brought to the laboratory. Samples of internal organs and yolk sacs are taken aseptically from the embryos and pooled.
	<ul style="list-style-type: none"> Chick Box Liners 	<ul style="list-style-type: none"> • Collect at least five basket liners per flock, total surface area at least 1 m², in a sterile bag and put this sample in a second bag before sending it to the laboratory.
Bacteria Isolation and Identification	<ul style="list-style-type: none"> Biological specimen from animals (Organ/tissue, faeces, fluid such as blood and urine, and other discharges) 	<ul style="list-style-type: none"> ➤ Fresh samples should be collected aseptically and suitable portions of each sample should be placed in an appropriate individual sterile container/petri dish, syringe, tube or bottle and clearly labelled. The specimen collected should preferably be chilled. Do not FREEZE. ➤ Swabs derived from the specimen are also acceptable. Collect using sterile swab tips and place in tube or bottle with BPW. Keep in a cool box.
Detection and Enumeration of Total Coliform / <i>Escherichia coli</i> (<i>E.coli</i>) in water	<ul style="list-style-type: none"> Drinking Water 	<ul style="list-style-type: none"> ➤ Collect at least 100 ml of water in sterile polypropylene sample containers with leak-proof lids (kept chilled during transport and storage). Use insulated containers to assure proper maintenance of storage temperature. Ensure that sample bottles are not totally immersed in water from melted ice during transit or storage.

Test to be performed	Appropriate samples	Collection of samples
Aerobic Plate Count for Environmental Monitoring (Air quality)	<ul style="list-style-type: none"> • Exposure plates 	<ul style="list-style-type: none"> ➤ Clearly mark the base of each plate date & time, area & location using permanent marker. Place the plates in some position with the lids still on. Remove the lid only once the sample is in place. Leave the plates exposed for 15 to 60 minutes. After exposing the plates replace carefully without touching the surface of the agar. Store the plates with the details in zip-lock bag include a control plate that has not been exposed.
Aerobic Plate Count for Environmental Monitoring (Surface)	Surface monitoring by swabbing/agar contact plate	<ul style="list-style-type: none"> ➤ Remove cover of the RODAC plate and gently push the agar surface on a level surface to examine for at least 10 seconds. Cover the plate with the lid. Store the plates with details in zip-lock bag.
Heterotrophic Plate Count for Environment Monitoring (Water)	Water	<ul style="list-style-type: none"> ➤ Collect at least 100 ml of water in sterile polypropylene sample containers with leak-proof lids (kept chilled during transport and storage). Use insulated containers to assure proper maintenance of storage temperature. Ensure that sample bottles are not totally immersed in water from melted ice during transit or storage.

The above transport media are available (by self-collection) from the VLS. Please contact VLS reception at least two days in advance.

NOTE: All samples for microbiology culture should be collected **PRIOR TO ANTIBIOTIC TREATMENT** (unless for monitoring treatment efficacy). It is highly recommended that collecting all these samples should be taken aseptically (wear gloves), placed them securely in sterile leak-proof containers / bags / bottles and transported in chilled condition directly to the VLS on the day of collection.

MOLECULAR BIOLOGY

1. Samples for Avian Influenza Confirmation Testing include:

- Swabs (nasal, tracheal, cloacal)
 - Whole blood and body fluids
-
- Samples should be kept chilled in a suitable cool box and submitted to the laboratory without any delay and clearly labeled with sample identification and date of collection.
 - In case of a swab, use a sterile recommended cotton swab, where the tip of it has to swab vigorously to exfoliate cells. Place the swab immediately after collecting into sterile tube without any additives provided by the laboratory.
 - Whole blood samples can be sent in syringes used for collecting or aliquoted into vials or tubes with each containing volume not less than 1.0 ml.

2. Animal raw, processed meat and animal feed samples for species identification / meat speciation testing.

- Raw, processed meat or animal feed samples should be collected at least 50 g each.
- Each sample should be placed in an individual sealed bag or container and labelled.

3. Samples for Rabies Testing include:

- Saliva
 - Oral swab
-
- a) For swab: Use sterile cotton swab to swab the anterior surface of tongue and check mucosa of dog for 20 s. Immersed the swab in PBS (2 mL) and kept in cool box filled with ice pack. Label tube correctly with sample identification and date of collection and send to VLS on the same day of sample collection. Swab and PBS are provided by VLS.

- b) For saliva: Using 10 or 50 ml sterile syringe, collect about 0.5 ml dog saliva and dispense it into sterile tube. Label tube correctly and send to VLS on the same day of sample collection.

CLINICAL CHEMISTRY / HAEMATOLOGY

Clinical Chemistry / Haematology laboratory currently carries out three types of tests: (i) for blood chemistry, (ii) complete blood cell count (CBC) and (iii) urine analysis (UA) of both large and small animals.

Blood collection tubes with K₂ EDTA anticoagulant are required to take samples for use in haematology tests, whereas those with lithium heparin ones are used to take blood for use in biochemistry examinations (blood chemistry tests). Both types of tubes are not vacutainers.

Whole blood samples should be collected using syringes, and collected blood samples are immediately transferred into respected tubes that are provided by the laboratory and should be labelled clearly with animal identification and date of collection. Each sample should be taken at least 1 ml in volume or equal to a half-full mark filled of a 3 ml similar kind tube. Please make sure that blood withdrawn in the tubes is mixed well immediately with the anticoagulant by a 10 time-gently inversion of the tubes.

Urine samples should be collected in new-clean urine tubes provided. All samples collected should be clearly and precisely labeled with animals' identification and date of collection by using a permanent marker and kept in a cool box with ice or frozen ice packs to control temperature not exceeding 8°C during submission to the laboratory on the same day of collection.

Urine samples may also be collected fresh and kept in new-clean urine tubes provided. The tube should also be clearly and precisely labelled. The sample can be submitted to the laboratory without keeping it in cool box, however it must be kept in clean seal plastic bag and submitted to the laboratory upon collection.

TABLE 1: GUIDELINES FOR TAKING SAMPLE SIZES OF ANIMAL POPULATIONS

Note: Guidelines shown below is for the Diagnostic's sample only.

Sample size required providing 95% certainty of detecting at least one positive case at a prevalence of 5%; and size required in a population where there is a known disease. (Adopted from Cannon and Roe, 1986)

Population Size (N)	Sample size set based on a 5% prevalence (n/N x 100%)	Sample size required where infectious disease shown on being suspected*
<10	All of the samples	4
10	10	4
20	19	4
30	26	4
40	31	5
50	35	5
60	38	5
70	40	5
80	42	5
90	43	5
100	45	5
120	47	5
140	48	5
160	49	5
180	50	5
200	51	5
250	53	5
300	54	5
350	54	5

* with a prevalence rate as at least 50%

Population Size (N)	Sample size set based on a 5% prevalence (n/N x 100%)	Sample size required where infectious disease shown on being suspected*
400	55	5
450	55	5
500	56	5
600	56	5
700	57	5
800	57	5
900	57	5
1000	57	5
1200	58	5
1400	58	5
1600	58	5
1800	58	5
2000	58	5
3000	58	5
4000	58	5
≥5000	59	5
∞	59	5

Note: n= Number of animals expected diseased or shown clinical signs

N= Total Population (population at risk)

SAMPLE RETENTION TIME

All samples received will be used for the whole process of testing. Any leftovers will be retained by respective laboratories at most 2 weeks after Laboratory Analysis Report is released.

CLIENT CHARTER

The table below indicates the turnaround time required for each tests.

Table 2: Client charter

No.	Laboratory	Laboratories activities/ Scope of work (Analyses)	Terms (in No. of working days)
1	Pathology	To Perform Pathological Examinations/ Diagnosis	3
2	Parasitology	To Perform Parasitological Examinations	3
3	Veterinary Microbiology	To Perform Veterinary Microbiological Testing	12
4		To Perform Environmental Monitoring Testing	4
5		To Perform Environmental Monitoring Testing (<i>Salmonella</i>)	12
6		To Perform Detection and Enumeration of Total Coliform and <i>E.coli</i> in Water	3
7	Serology	To Perform Serological Test using Rapid Test	1
8		To Perform Serological Test using ELISA Test	5
9	Molecular Biology	To Perform Molecular Biological Testing	5
10	Clinical Chemistry/ Haematology	To Perform Clinical Chemistry Testing	3

The VLS shall always try to produce reports as soon as any analyses are finished.

Every customer shall collect his/her reports of analysis from the VLS counter.

All valued customers are advised to fully understand on terms of every client charter for their required analyses, this is to avoid any misunderstanding on analysis report collection from the VLS counter desk. Your cooperation is greatly appreciated.

NEW CLIENTS

New customer is required to fill up the New Testing Request Form, VLS-SOP-09-02 and attend a short briefing on the tests offered by VLS. Existing customers requesting for new test are also required to fill the above form. The form is available upon request at the reception counter.

**FOR FURTHER ENQUIRIES
PLEASE CALL / CONTACT
VETERINARY LABORATORY SERVICES UNIT**

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